

## **GUIDE TO THE ANALECT SPECTRAL FILE FORMAT**



Applied Instrument Technologies, Inc.  
2121 Aviation Dr  
Upland, CA 91786  
Phone: (909)204-3700



Guide to the Analect Spectral File Format

DOC. NO.: 98-0451

**DOCUMENT RELEASE RECORD (DRR)**

RELEASES (& DATES)	Update this document by:			AUTHORIZATIONS (& REMARKS)
	ACTIVE PAGES	PAGES ADDED	PAGES DELETED	
Issue-1 (4/20/2004)	All.	All.	None.	Initial release.
Issue-2 (5/4/2006)	All.	ii	None.	SPR#1280 Added document release record as page ii. Updated document information and indices throughout. Corrected document number in Appendix and Index. Updated division name on cover page and in section 1.1. Added RPM product line reference in section 1.2. Corrected figure 1 alignment and caption in section 2. Corrected abbreviation definition on first use throughout. Added new section 5.
Issue-3 (6/20/2006)	All.	None.	None.	SPR#1290 Corrected headers in Appendix & Index, updated field codes and issue number.
Issue-4 (9/28/2012)	All	None	None	SPR #122 Company Name and Logo updated on front cover page.



## TABLE OF CONTENTS

<b>1</b>	<b>INTRODUCTION.....</b>	<b>1</b>
1.1	Purpose .....	1
1.2	Scope .....	1
1.3	Definitions, Acronyms, and Abbreviations .....	1
<b>2</b>	<b>BASIC STRUCTURE.....</b>	<b>2</b>
<b>3</b>	<b>DESCRIPTORS .....</b>	<b>3</b>
<b>4</b>	<b>ASF HEADERS.....</b>	<b>4</b>
<b>5</b>	<b>HEADER VERSIONS AND RAMAN REDEFINITIONS .....</b>	<b>5</b>
<b>6</b>	<b>VARIOUS TYPE DEFINITIONS .....</b>	<b>7</b>
	<b>APPENDIXES.....</b>	<b>A-1</b>
	Appendix A - Definitions, Acronyms, and Abbreviations .....	A-1
	<b>INDEX.....</b>	<b>I-1</b>

## LIST OF FIGURES

Figure 1 - ASF file block diagram .....	2
---	---

## LIST OF TABLES

Table 1 - Raman redefinitions of ASF header parameters.....	5
---	---



---

## **1 Introduction**

### **1.1 Purpose**

The purpose of this document is to serve as an overview of the Analect Spectral File (ASF) format used by the Applied Instrument Technologies (AIT), Upland, California. This document is identified as the Guide to the Analect Spectral File Format, document number 98-0451.

### **1.2 Scope**

The Analect spectral file (ASF) format is the primary data file format used for the Analect and RPM product lines. Providing a clear and concise description of the file format allows for quicker and more reliable internal and external software development. This guide describes the overall structure of the Analect spectral file format, and details the individual elements and definitions used within the file header and various descriptors.

### **1.3 Definitions, Acronyms, and Abbreviations**

Please refer to Appendix A for clarification of terminology used in this document.



## 2 Basic Structure

The basic structure of an Analect spectral file consists of a number of descriptor blocks that indicate the type, size, and location of a corresponding file component. In addition, each descriptor block in a file contains the location of the next descriptor block and thus the next component. The following diagram illustrates the overall structure of an Analect spectral file.




Figure 1 – Analect spectral file block diagram



### 3 Descriptors

A descriptor is a 16 byte block that indicates a components type, size, and location in the file. The structure of a descriptor block is as follows:

```
typedef struct
{
    long ld;    // link down, position of next descriptor
    long la;    // link across, unused at this time
    long size;  // size of corresponding component, including descriptor
    int version; // version of ASF file structure
    char ctype; // component type
    char ftype; // file type
} DESCRIPTOR;
```

The component type is defined as follows:

```
#define UNDEFINED 0
#define TRACEDATA 1
#define TRACEHEADER 2
#define PEAKTABLE 3
#define COMMENT 4
#define CMDHISTORY 5
#define AFHEADER 6
```

The file type is defined as follows:

```
#define UNDEFINED 0
#define TRACEFILE 1
#define GCFILE 2
#define PARVFILE 3
#define DATABASEFILE 4
```

## 4 ASF Headers

An ASF header is an 898 byte block containing information about the collection of the ASF data that follows later in the file. The structure of the ASF header is as follows:

```
typedef struct
{
    long
        time,                // time elapsed in seconds since 00:00:00 Jan 01 1970.
        serial_no,           // serial number
        ndata,               // # of trace data points
        ig_size,             // # of points in original interferogram
        fft_size,            // # of points in interferogram fast Fourier transform (FFT)
        fft_spin,            // # of points skipped at beginning of transform result
        scans_sig,           // # of signal scans
        scans_bkg,           // # of background scans
        _L_FI_LL[ 6];        // room for additional longs in the future
    float
        xleft,              // x-axis leftmost value
        xright,             // x-axis rightmost value
        yorg,               // y-axis default origin for display (center of full scale)
        ymax,               // y-axis default maximum for display
        yscale,             // y-axis scale factor for integer format trace functions
        ig_step,            // retardation step for points in original interferogram
        resolution,         // spectral resolution in wave numbers(WN)
        mol_wt,             // molecular weight in atomic units (A. U.)
        bp,                 // boiling point (degrees Celsius)
        mp,                 // melting point (degrees Celsius)
        xdelta,             // x-axis delta between points, = (xright-xleft)/(ndata-1)
        laser_wn,           // laser wave number for Raman instruments
        _F_FI_LL[ 2];       // room for additional floats in the future
    int
        lgain_sig,          // 2*log2(signal gain)
        lgain_bkg,          // 2*log2(background gain)
        phi_glen,           // size of phase interferogram if Mertz method
                           // phase-corrected spectrum or Foreman method
                           // phase-corrected interferogram
        ver_num,            // header structure version number
        transept,           // flag word for transept interferometer data,
                           // bit 0 set if non-linear spectrum or igram
        pc_flags,           // phase correct flag word, see definitions
        _I_FI_LL[ 6];       // room for additional integers in the future
    TFMF_TYPE trace_fmt;    // trace format, see definitions
    DFMT_TYPE data_fmt;     // trace data format, see definitions
    XAXIS_TYPE xaxis;       // x-axis units, see definitions
    YAXIS_TYPE yaxis;       // y-axis units, see definitions
    BSPLT_TYPE bs_type;     // beamsplitter and wedge type, see definitions
    APOD_TYPE ap_type;      // apodization type, see definitions
    int _E_FI_LL[ 2];       // room for additional enumerations in the future
    char
        title[ 60],         // trace title for display and plot
        desc1 [ 60],        // trace description lines, usually brief
        desc2 [ 60],        // sample description if spectrum trace
        mfg [ 24],          // instrument manufacturer
        model [ 24],         // instrument model
        origin [ 60],       // trace origin
        owner [ 60],        // trace (copyright) owner
        operator [ 60],     // instrument operator
        casnumber [ 16],    // Chemical Abstracts Service (CAS) registry number
        casname [ 60],      // CAS name
        mol_form [ 60],     // molecular formula
        wws [ 32],          // Wsweiser index
        xunits [ 8],        // x-axis display units

```

```

yunits [ 8],          // y-axis display units
detector [ 16],       // detector type
int_type [ 16],       // interferometer type
ap_comm [ 26],        // apodization type annotation
_C_FILL[ 96];         // room for additional characters in the future
} ASFHEADER;

```

## 5 Header Versions and Raman Redefinitions

In order to maintain compatibility with the existing header structure, several of the data items defined in the ASFHEADER structure shown in section 4 are redefined for Raman systems.

The ASF header contains an internal header structure version number in the integer “ver\_num”. The current header structure version is 3.10, and is defined in “ver\_num” as the integer value 310. Any spectra with header version numbers less than 310 are considered to contain Fourier transform infrared (FTIR) data.

Any spectra with header version numbers greater than or equal to 310 may contain either FTIR or Raman data. The floating point value “laserwn” is new with version 3.10, and is used to determine if the data is from an FTIR or Raman system. Header structures prior to version 3.10 were the same size, but the number of spare floats was 3 instead of 2, and there was no “laserwn” floating point value.

If the “ver\_num” is greater than or equal to 310 and the “laserwn” value is between 9,400 and 50,000 inclusive, then the data is considered to be Raman. Otherwise, the data is considered to be FTIR.

Since the display of FTIR data items follows the header structure parameters, the items redefined for Raman systems are shown in the following table.

Table 1 - Raman redefinitions of ASF header parameters

ASF Header Parameter	Raman Meaning
Title	Acquisition information
Desc1	Comment
Desc2	X correction information
ScansSig	Exposures co-added
Wws	Exposure period, milliseconds (ms)
Igstep	ASF point spacing, inverse centimeters (cm <sup>-1</sup> )
Fftsize	Grating period, line pairs (lp)/millimeter (mm)
Molwt	Grating blaze, nanometers (nm)
Mp	Camera Temp, degrees Celsius
Bp	Camera Temp locked 0 or 1
Inttype	Spectrograph serial number
Laserwn (was _F_FILL[0])	Laser wavenumber



The “title” entry in Raman spectral file headers, as shown in Table 1, is filled with acquisition information. The acquisition information is further broken out as follows:

“S=” The number of the strip that produced the data.

“AQ=” The name of the acquire parameter file that directed the acquisition and calibration of the data.

“F=” A nine character code detailing the corrections performed upon the data defined as follows:

1st: Dark correct: N = none, F = File based, A = Automatic.

2nd: X correction performed on this spectrum true or false (T or F).

3rd: X correction data obtained from this spectrum T or F.

4th: Y correction performed T or F.

5th-9th: X correction points N or 0 or 1. For this spectrum the neon correction list is examined if a point is defined then a 0 or 1 is entered. If no point is defined, a N is entered.

“%F=” The maximum signal in percent of full scale.

An example of a Raman title entry is:

S=3 AQ=N1S\_30Z F=FTTT11111 %F=24.2%

The “desc2” entry in Raman spectral file headers, as shown in Table 1, is filled with X correction information. The X correction information is further broken out as follows:

“RA=” Average offset between known and measured Raman references.

“LO=” Laser offset correction from AQP file.

“A0=” Constant term of quadratic correction equation.

“A1=” Proportional term of quadratic correction equation.

“A2=” Square term of quadratic correction equation.

## 6 Various Type Definitions

The trace format is defined as follows:

```
typedef enum TFMT_TAG
{
    TFMT_UNK,      // unknown or undefined
    TFMT_SPC,      // signal spectrum (ratioed)
    TFMT_BKG,      // background or single-beam spectrum
    TFMT_IGM,      // interferogram trace
    TFMT_TLM,      // time function trace
    TFMT_CGM,      // chromatogram
                  // gas chromatogram (GC)
                  // routine gas chromatogram (RGC)
                  // liquid chromatogram (LC), etc.
    TFMT_ARB,      // arbitrary trace function
} TFMT_TYPE;
```

The trace data format is defined as follows:

```
typedef enum DFMT_TAG
{
    DFMT_UNK,      // unknown or undefined (no data)
    DFMT_INT2,     // 2-byte (16-bit) signed binary integer data
    DFMT_INT4,     // 4-byte (32-bit) signed binary integer data
    DFMT_INT8,     // 8-byte (64-bit) signed binary integer data
    DFMT_FLT4,     // 4-byte (32-bit) IEEE floating point data
    DFMT_FLT8,     // 8-byte (64-bit) IEEE floating point data
} DFMT_TYPE;
```

The x-axis units is defined as follows:

```
typedef enum XAXIS_TAG
{
    XAX_UNK,      // unknown or undefined x units
    XAX_WN,       // trace function of frequency, normally in wavenumbers (cm-1)
    XAX_MCR,      // trace function of wavelength, normally in microns
    XAX_TIME,     // trace function of time
    XAX_ARB,      // trace function of arbitrary x units
} XAXIS_TYPE;
```

The y-axis units is defined as follows:

```
typedef enum YAXIS_TAG
{
    YAX_UNK,      // unknown or undefined y units
    YAX_TR,       // transmittance
    YAX_AB,       // absorbance (or kumu)
    YAX_PAS,      // photoacoustic trace
    YAX_ARB,      // arbitrary y-axis units
} YAXIS_TYPE;
```

The beamsplitter and transept wedge composition are defined as follows:

```
typedef enum BSPLT_TAG
{
    BSPLT_UNK,    // unknown beamsplitter type
    BSPLT_KBR,    // potassium bromide (KBr)
    BSPLT_CSI,    // cesium iodide (CsI)
    BSPLT_CAF,    // calcium fluoride/fluorite (CaF2)
    BSPLT_MYLAR, // mylar [(r) DuPont]
```



Guide to the Analect Spectral File Format

---

```
BSPLT_NACL, // sodium chloride (NaCl)
BSPLT_BAF  // barium fluoride (BaF2)
} BSPLT_TYPE;
```



## Guide to the Analect Spectral File Format

The apodization is defined as follows:

```
typedef enum APOD_TAG
{
    APDT_UNK,    // unknown or not applicable
    APDT_BOX,    // Boxcar (no apodization)
    APDT_NBWEAK, // Norton-Beer (weak)
    APDT_NBMED,  // Norton-Beer (medium)
    APDT_NBSTRT, // Norton-Beer (strong)
    APDT_GAUSS,  // Gaussian
    APDT_HG,     // Happ-Genzel
    APDT_TRAP,   // Trapezoidal
    APDT_TRI,    // Triangular
    APDT_TRI2,   // Triangular squared
    APDT_BESS,   // Bessel
    APDT_COS,    // Cosine
    APDT_SINC2,  // Sinc squared
    APDT_BH3T,   // 3 term Blackmann-Harris
    APDT_BH4T    // 4 term Blackmann-Harris
} APOD_TYPE;
```

The bit definitions for the transept flag word are as follows:

```
#define TSI_NONLIN 0x0001 // File is nonlinear (interferogram or spectrum)
#define TSI_TSI    0x0002 // File taken on an Analect Transept class
                        // interferometer.
```

The bit definitions for the phase correct flag word are as follows:

```
#define PCF_MSHFT 0 // Bits 0-3 are method
#define PCF_MMASK 0x000F
#define PCM_SSMERTZ 0 // 0 = Single-sided Mertz Method
#define PCM_SSMAG 1 // 1 = Single-sided Magnitude Method
#define PCM_DSMAG 2 // 2 = Double-sided Magnitude Method
#define PCM_SSFOR 3 // 3 = Single-sided Forman Method
#define PCF_TSHFT 4 // Bits 4-7 are truncation type
#define PCF_TMASK 0x00F0 // (Reserved for future implementation)
#define PCT_RAMP 0 // 0 = Ramp, over width of phase igram
#define PCT_RECT 1 // 1 = Rectangular truncation at
                        // zero path difference (ZPD)
```

**Appendixes****Appendix A - Definitions, Acronyms, and Abbreviations**

A.U.	Atomic units.
AIT	Applied Instrument Technologies.
ASF	Analect spectral file.
BaF2	Barium fluoride.
CaF2	Calcium fluoride.
CAS	Chemical abstracts service.
cm-1	Inverse centimeters.
CsI	Cesium iodide.
FFT	Fast Fourier transform.
FTIR	Fourier transform infrared.
GC	Gas chromatogram.
SLS	Space, Land & Sea.
KBr	Potassium bromide.
LC	Liquid chromatogram.
lp	Line pairs.
ms	Milliseconds.
NaCl	Sodium chloride.
nm	Nanometers.
RGC	Routine gas chromatogram.
WN	Wavenumber.
ZPD	Zero path difference.

## Index

---

### *A*

ASF Headers · 4

---

### *B*

Basic Structure · 2

---

### *D*

Definitions, Acronyms, and Abbreviations ·  
A-1

Descriptors · 3

---

### *H*

Header Versions and Raman Redefinitions ·  
5

---

### *I*

Introduction · 1

---

### *P*

Purpose · 1

---

### *S*

Scope · 1

---

### *V*

Various Type Definitions · 7